



solutions for welfare





We specialize in wholesale of products that enable improved quality of life and daily activities of elderly, rugged, disabled and health professionals, offering a trusted orthopedic products and high quality.

Our products are approved by the **Ministry of Health and Consumer** licensed manufacturer of medical devices N° 4643-PS.

Products that meet the requirements set by the European Economic Community on health products Class I.

Wide range of geriatrics, orthopedics and pharmacy products.

Spanish manufacture products that stand out for their very competitive design and quality-value ratio.



676 767 176




dioni@comercialnespral.com



comercialnespral.com





Crutches

4



Collars

5



Slings

6



Restraining Belts

7



Immobilizers

8



Cushions

10



Protectors

11



Anatomical Crutch

D-01030

Anatomical Crutch extendible lower part with safety clip.

Made of aluminum tube and adjustable in height by metal stent. The handles and elbow pads are made of polypropylene. The ferrule is P.V.C. with draining system and with the same grip.

Max weight: 125 kg
Max arm length: 225 mm
Min arm length: 225 mm
Handle length: 125 mm
Support angle: 21°
Grip angle: 104,5°
Cane length: 674-980 mm

Anatomical Crutch

D-01050

Anatomical Crutch extendible lower part with safety clip.

Made of aluminum tube and adjustable in height by metal stent. The handles and elbow pads are made of polypropylene. The ferrule is P.V.C. with draining system and with the same grip.

Max weight: 125 kg
Max arm length: 230 mm
Min arm length: 230 mm
Handle length: 125 mm
Support angle: 21°
Grip angle: 104,5°
Cane length: 630-980 mm



Anatomical Extendible Crutch

D-01070

Anatomical Crutch extendible in upper and lower part with safety clip.

Made of aluminum tube and adjustable in height by metal stent. The handles and elbow pads are made of polypropylene. The ferrule is P.V.C. with draining system and with the same grip.

Max weight: 125 kg
Max arm length: 250 mm
Min arm length: 200 mm
Handle length: 125 mm
Support angle: 21°
Support angle: 104,5°
Cane length: 674-980 mm

Articulate Extendible Crutch

D-01090

Articulate crutch extendible in upper and lower part with safety clip.

Made of aluminum tube and adjustable in height by metal stent. The handles and elbow pads are made of polypropylene. The ferrule is P.V.C. with draining system and with the same grip.

Max weight: 125 kg
Max arm length: 250 mm
Min arm length: 200 mm
Handle length: 125 mm
Support angle: 21°
Support angle: 104,5°
Cane length: 674-980 mm






Polyethylene Collar

D-03010

To assure lightness and strength, perforated in order to guarantee comfort and hygiene. The lower and upper edges are lined with soft rubber and artificial leather, fastened with elastic tape.

INDICATIONS: Light or moderated cervical sprains without stability, soft areas traumatism and other light cervical injuries, rheumatic processes (inflammatory or degenerative).

COLOUR: 

HEIGHT: 8-10 cm (adjustable).

SIZE:

S: 44 cm length.

M: 48 cm length.

L: 53 cm length.

XL: 57 cm length.

XXL: 62 cm length.

Foam Collar

D-03030

Collar whose base is a polyurethane foam D30 inside and outer tubular cotton, fastened with elastic tape.

Anatomically shaped.

INDICATIONS: Light or moderated cervical sprains without stability, soft areas traumatism and other light cervical injuries, rheumatic processes (inflammatory or degenerative).

COLOUR: 

HEIGHT: 9 cm.

SIZE:

S: 40 cm length.

M: 44 cm length.

L: 49 cm length.

XL: 54 cm length.

XXL: 59 cm length.



Semi-Rigid Foam Collar

D-03035

Collar whose base is a polyethylene foam uncrosslinked tubular inside and outside cotton, fastened with elastic tape.

Anatomically shaped.

INDICATIONS: Light or moderated cervical sprains without stability, soft areas traumatism and other light cervical injuries, rheumatic processes (inflammatory or degenerative).

COLOUR: 

HEIGHT: 9 cm.

SIZE:

S: 40 cm length.

M: 44 cm length.

L: 49 cm length.

XL: 54 cm length.

XXL: 59 cm length.

Mentalis Support Collar

D-03040

Collar whose base is a polyethylene hemming and padded with soft rubber and artificial leather, fastened with elastic tape with support for chin on top.

INDICATIONS: Light or moderated cervical sprains without stability, soft areas traumatism and other light cervical injuries, rheumatic processes (inflammatory or degenerative).

COLOUR: 

HEIGHT: 8-10 cm (regulable).

SIZE:

S: 44 cm length.

M: 48 cm length.

L: 53 cm length.

XL: 62 cm length.





Closed Sling

D-02140

INDICATIONS: Download the shoulder and collarbone, inflammatory processes which require restraint, clavicle and humerus fractures, holding the arms in casts, sprains and dislocations of the upper limb.

HEIGHT: Adjustable | **SIZE:** S / M / L / XL | **COLOURS:** ☐ ☒

Sling Tape

D-02142

INDICATIONS: Download the shoulder and collarbone, inflammatory processes which require restraint, clavicle and humerus fractures, holding the arms in casts, sprains and dislocations of the upper limb.

HEIGHT: Adjustable | **Only one size** (115x40 mm) | **COLOURS:** ☐ ☒



Shoulder Immobilizer

D-02145

Closed sling with shoulder immobilizer, made of padded cloth allowing for the adherence of Velcro on the entire surface. It may be adapted to patients of all ages, and may be used for both left and right sides.

INDICATIONS: Download the shoulder and collarbone, inflammatory processes which require restraint, clavicle and humerus fractures, holding the arms in casts, sprains and dislocations of the upper limb.

Only one size (200x40x350 mm) | **COLOURS:** ☒ ☐

Sling in Ratier Tape

D-02146

Allows front or back posture and under the arm.

INDICATIONS: Download the shoulder and collarbone, inflammatory processes which require restraint, clavicle and humerus fractures, holding the arms in casts, sprains and dislocations of the upper limb.

HEIGHT: Adjustable | **Only one size** | **COLOURS:** ☐ ☒





Seat Restraining Belt

INDICATIONS: The main therapeutic target of this product range is to provide good stability and user safety wheelchair or geriatric chair so as to prevent the patient can slip or move sideways voluntarily. At the same time, be maintained in the patient's spine remain erect position. They are indicated in very hypotonic patients. Not completely immobilizes the user but allows you to make small movements around the equilibrium position. Depending on the models, the product works on a bigger or smaller body surface.

Only one size | COLOUR: ●

D-2100

Seat with Pelvic Restraining Belt

INDICATIONS: The main therapeutic target of this product range is to provide good stability and user safety wheelchair or geriatric chair so as to prevent the patient can slip or move sideways voluntarily. At the same time, be maintained in the patient's spine remain erect position. They are indicated in very hypotonic patients. Not completely immobilizes the user but allows you to make small movements around the equilibrium position. Depending on the models, the product works on a bigger or smaller body surface.

D-2101

Only one size | COLOUR: ●



Strapped Seat Restraining Belt

INDICATIONS: The main therapeutic target of this product range is to provide good stability and user safety wheelchair or geriatric chair so as to prevent the patient can slip or move sideways voluntarily. At the same time, be maintained in the patient's spine remain erect position. They are indicated in very hypotonic patients. Not completely immobilizes the user but allows you to make small movements around the equilibrium position. Depending on the models, the product works on a bigger or smaller body surface.

Only one size | COLOUR: ●

D-2102


Strapped Seat with Pelvic Restraining Belt

INDICATIONS: The main therapeutic target of this product range is to provide good stability and user safety wheelchair or geriatric chair so as to prevent the patient can slip or move sideways voluntarily. At the same time, be maintained in the patient's spine remain erect position. They are indicated in very hypotonic patients. Not completely immobilizes the user but allows you to make small movements around the equilibrium position. Depending on the models, the product works on a bigger or smaller body surface.

D-2103

Only one size | COLOUR: ●





Abdominal Magnet Immobilizer

THERAPEUTIC PURPOSES: The aim of this product family is to immobilize or limit the patient's movements, avoiding dangerous falls or moves from the bed. Depending on the model and accessories used certain postural changes are possible. They are suitable especially for the immobilization of people with impaired mental faculties, to prevent falls, catheter, etc.

Long (including strips): 260 mm. Length: 135 mm. Width: 18 mm

Only one size | COLOUR: ☐

D-2040

Extremities magnet immobilizer

THERAPEUTIC PURPOSES: The aim of this product family is to immobilize or limit the patient's movements, avoiding dangerous falls or moves from the bed. Depending on the model and accessories used certain postural changes are possible. They are suitable especially for the immobilization of people with impaired mental faculties, to prevent falls, catheter, etc

D-2050

Only one size | COLOUR: ☐



Extremities Magnet Immobilizer

THERAPEUTIC PURPOSES: The aim of this product family is to immobilize or limit the patient's movements, avoiding dangerous falls or moves from the bed. Depending on the model and accessories used certain postural changes are possible. They are suitable especially for the immobilization of people with impaired mental faculties, to prevent falls, catheter, etc.

Only one size | COLOUR: ☐

D-2060

Thorax Magnet Immobilizer

THERAPEUTIC PURPOSES: The aim of this product family is to immobilize or limit the patient's movements, avoiding dangerous falls or moves from the bed. Depending on the model and accessories used certain postural changes are possible. They are suitable especially for the immobilization of people with impaired mental faculties, to prevent falls, catheter, etc.

D-2090

Only one size | COLOUR: ☐





Abdominal Immobilizer Standard

THERAPEUTIC PURPOSES: The aim of this product family is to immobilize or limit the patient's movements, avoiding dangerous falls or moves from the bed. Depending on the model and accessories used certain postural changes are possible. They are suitable especially for the immobilization of people with impaired mental faculties, to prevent falls, catheter, etc.

Only one size | COLOUR: ☐

D-2010

Extremities Immobilizer Standard

THERAPEUTIC PURPOSES: The aim of this product family is to immobilize or limit the patient's movements, avoiding dangerous falls or moves from the bed. Depending on the model and accessories used certain postural changes are possible. They are suitable especially for the immobilization of people with impaired mental faculties, to prevent falls, catheter, etc.

D-2020

Only one size | COLOUR: ☐



Universal Knee Splints

THERAPEUTIC PURPOSES: The aim of this product family is to immobilize or limit the patient's movements, avoiding dangerous falls or moves from the bed.

SIZES	Body	S 600x400x6	M 500x350x6	L 400x300x6
(mm)	Plate	S 540x20x3	M 400x20x3	L 360x20x3

In 3 adjustable panels | COLOUR: ☒

D-2500

Clavicle Immobilizer

THERAPEUTIC PURPOSES: The aim of this product family is to immobilize or limit the patient's movements, avoiding dangerous falls or moves from the bed. Depending on the model and accessories used certain postural changes are possible. They are suitable especially for the immobilization of people with impaired mental faculties, to prevent falls, catheter, etc.

D-2600

SIZES: S / M / L | COLOUR: ☒





Anti-Eschar Cushion

D-6300

SQUARE CUSHION.

INDICATIONS: Prevention and treatment of low back pain and degenerative processes. Occupational pathology. Maintenance of body heat, etc.

Only one size | COLOUR: ☐

Anti-Eschar Cushion

D-6350

SQUARE CUSHION WITH HOLE.

INDICATIONS: Prevention and treatment of low back pain and degenerative processes. Occupational pathology. Maintenance of body heat, etc.

Only one size | COLOUR: ☐



Anti-Eschar Cushion

D-6400

ROUND CUSHION WITH HOLE.

INDICATIONS: Prevention and treatment of low back pain and degenerative processes. Occupational pathology. Maintenance of body heat, etc.

Only one size | COLOUR: ☐

Anti-Eschar Cushion

D-6450

HORSESHOE CUSHION.

INDICATIONS: Prevention and treatment of low back pain and degenerative processes. Occupational pathology. Maintenance of body heat, etc.

Only one size | COLOUR: ☐



Cushions



Viscoelastic Cushion

D-6200

Memory foam cushion D50 in a single block. 3D fabric cover with non-slip.

INDICATIONS: Prevention and treatment of low back pain and degenerative processes. Occupational pathology. Maintenance of body heat, etc.

Only one size | COLOUR: ☒

Protectors



Anti-Eschar Protector

D-2700

FOOT OR ELBOW PROTECTOR DELUXE.

Only one size (100x100x100 mm) | COLOUR: ☐

Anti-Eschar Protector

D-2701

FOOT OR ELBOW PROTECTOR STANDARD.

Only one size (100x100x100 mm) | COLOUR: ☐



